

SHERA  
PURE

Super hard stone



easily readable · colour reduced  
non - splintering · extreme hardness

First it is creamy, later it is extremely hard – the dental stone SHERAPURE offers best technical properties in pure form. It has a maximum of dimensional stability, it is absolutely precise in reproducing details and very harmonic thanks to the discreet colors. SHERAPURE allows your work to leave an impressive effect – very pure created by Andreas Nolte.

Created by  Andreas Nolte

# SHERA PURE

## Super hard stone

Visionary in dental technology Andreas Nolte has shaped SHERAPURE and has incredibly supported the development. Get to know our latest dental stone and appreciate it as we do!

For implant work, sectioned, master,  
static and control models

Available in



Item no.	Unit
101155-25 (light grey) 101155-82 (cream white)	18 kg (6 x 3 kg)
101153-25 (light grey) 101153-82 (cream white)	3 kg sachet

### Technical specifications

mixing ratio (p : dw)*	100 : 20
hardness after 2 h (MPa; N/mm <sup>2</sup> )	310
hardness after 24 h (MPa; N/mm <sup>2</sup> )	330
compressive strength after 1 h (MPa; N/mm <sup>2</sup> )	70
compressive strength after 1 h (psi)	10,150
setting expansion after 2 h (%)	0.08
working time (min)	6.0 - 7.0
end of setting (min)	approx. 30